


ogawa
by OGAWA

(((pulse))) *Lite*



User Manual 使用手册

Model: OE 1031

CONTENTS

Safety Instructions	2
Warranty Claims	4
Specification	4
Name and Function of the Product	5
Program Mode	7
Use method	9
Storage and Maintenance	10
Surface cleaning	
Common Problems & Solutions	11
Contact us	12

Safety Instructions

Please read the manual carefully before using the appliance.
Power source for this appliance is 220–240V~ ;
alternating current frequency is 50/60Hz.

i. Usage environment

- To avoid malfunction or electric shock, do not use the appliance in a wet condition or dusty environment.
- Do not use the appliance in the bathroom.
- Do not use the appliance on the bed.
- Place the appliance in a safe location to avoid water immersion.
- Keep the appliance away from fire and heat, do not expose it to direct sunlight.
- Please use the appliance on a flat floor.

ii. Health precautions

If you are experiencing health complications or conditions as stated below, please consult your doctor before use:

- . Suffering from heart disease, malignancy, scoliosis, hypertension, osteoporosis, etc.;
- . Pregnancy or menstruation;
- . Fever, swelling, irritation, inflammation, congestion, allergy or any other health concerns

iii. When the product should not be used

- If the use exceeds 20 minutes at one time.
- This unit should not be used by children or invalids without adult supervision.
- Do not use this product when you are sleepy or drunk.
- If water is accidentally spilled onto the appliance.
- If the user feels abnormal or severe pain during the massage.
- If the user detects any malfunction or other abnormal conditions during operation.
- If there is lightning.

Safety Instructions

iv. Caution and warning

- Please completely insert the plug into outlet before use.
- Unplug the appliance from the outlet when not in use.
- Do not use any attachments that are not recommended by OGAWA.
- Do not operate this appliance if it has a damaged cord or plug. Please return the unit to OGAWA Service Center for examination and repairs.
- Do not drop or insert any object into any opening.
- Do not lift or carry this appliance by the power cord or any cord.
- Do not cover the appliance during operation.
- Do not crush the appliance and avoid any sharp folds.
- Anyone whose weight is more than 100KG should not use the appliance.
- This appliance is designed only for personal use, and not for professional purposes.
- The appliance can only be used by one person at a time.

v. Note

- The appliance is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge. They need to be supervised and instructed about its usage by a person responsible for their safety.
- Extensive usage could lead to excessive heating and a shorter product lifespan. Should this occur, discontinue the use and allow the unit to cool before operating again.
- If you are under medication or with a medical condition, please consult your doctor before use. Please do not use the product if you are not feeling well.
- Do not remove the cover of the appliance. Mishandling the internal components may result in malfunction or electric shock. Do not try to repair this appliance. The maintenance and repairs should only be made by technician.
- Do not connect or disconnect the power plug from the electrical outlet with wet hands.
- Do not use the appliance if the electrical outlet is loose.


Warranty Claims

1. Warranty is not transferable.
2. This warranty is effective only if the product is purchased from OGAWA or any authorized distributors.
3. Proof of purchase (original receipt) is required for all warranty repairs.
4. All implied warranties, including but not limited to these implied warranties of fitness and merchantability, are limited to the duration stated in the warranty policy.
5. Damage due to misuse, improper treatment and unauthorized modification and repair are not covered by this warranty.
6. Warranty does not cover accessories and attachments which do not belong to the product.
7. Warranty is not effective to rental, business, commercial, institutional or other nonresidential users.
8. All services covered by this warranty must be approved by OGAWA and repaired by authorized technicians only.
9. If replacement parts for defective materials are not available, OGAWA reserves the right to make substitution in lieu of repair or replacement.


Specification

Model No.	OE 1031
Power source	220~240V~
Alternating current	50/60Hz
Power consumption	200W
Max Timer	10 minutes
Weight	6.5KG
Dimension(L*W*H)	367x282x127mm

Name and Function of the Product



Name and Function of the Product



Power: Press the button with plug in, LCD shows "00" "10", the machine is switch on and standby.

Mode: Default to Auto program mode.

Auto program: "P1" "P2" "P3"

Start/Stop: Insert the power plug, switch on power, Then press the button to start the product. In the Operation of the machine, press this button to stop Working.

Speed+: Speed increase: Manual 1–30. Press this button. To increase the speed when the machine is operating.

Time+: Time increase: the default working time is 10 Minutes, divided into 10 levels, press this button to Increase time.



Speed-: Speed decrease: Manual 1–30. Press this button. To decrease the speed when the machine is operating.

Time-: Time decrease: the default working time is 10 Minutes, divided into 10 levels, press this button to Decrease time.




Note: The effective reception range of the remote Control signal is 2.5 meters, the infrared emission Indicator of the controller should be aligned the infrared Receiving window on the product.

Program Mode



Program Mode



Use method



Stand tiptoe in the middle of the massaging device. Straighten your back. Your abdominal muscles should be tight and you should feel changes in the tension in your calf muscles. You can also bend your knees 90 degrees. This stance will relax your waist and hips.

Standing on the massage device will align your feet and shoulders. Keep your back tight, your knees bent, and slowly tighten your leg muscles. To increase exercise intensity you should bend your knees as much as possible. You will begin to feel tension in your quads, hips, and back. This stance can exercise your thigh muscles, alleviate pain and soreness, promote blood circulation, and reduce fat deposits.



Kneel on the ground, place your hands on the product at shoulder width, and keep your body leaning forward. This stance will exercise your shoulder joints and arms.

Sit in front of the massager, with your back straight and your hands placed on either side. Place your calves on the massager, with your toes facing up. This stance should be used to completely relax while enjoying a powerful calf massage.



Sit on a chair in front of the massager. Place your feet on the device. The product's vibrations will help you relax the muscles in your back, buttocks, and thighs.

Storage and Maintenance

Surface cleaning

Before cleaning, ensure the power switch is turned off, and the attaching plug is pulled out of the socket.

The rubber pedal and covers can be cleaned with a wet cloth containing detergent, and the detergent should be removed with a Cloth



Notice

Do not clean the Product with corrosive detergent such as benzene Or thinner, as it may cause fading, corrosion or crack on the Product. Do not splash water into the Product to avoid short circuit.

Controller

As wet cloth may cause malfunction of the controller, only dry cloth Can be used.

Storage an Maintenance

Make sure that the power switch is turned off and the plug is removed From the power supply

Clean off the dust and dirt on the Product and keep it clean and away From places with high temperature or high humidity.

If the Product is not in use for along period, a dust proof cover (buy it Yourself) should be applied on the Product to prevent the entry of dust.

Common Problems & Solutions

S/N	Problems	Causes	Solutions
1	The Product can not be started	The plug is not in place	Insert the plug
		The power switch is off	Turn switch on, the indicator
		Th fuse is burn	Replace with a new fuse of the same mode
2	Sound like phew or Da-da is heard during working	Normal sound generated by the Motor or from the abrasion Between mechanical units	No additional treatment is required
3	Operate the controller, but the Product no action	The electricity of batteries is worn out	Replace batteries
		The controller is too far away from the Product	The effective distance of infrared signal is 2.5M
4	Sudden stop	Default working time is up	Restart it after 10 minutes' interval
5	Can not operation	Too many commands sent to the Product and they can not be responded in time	Power off and then restart

Note: Please consult after-sales department, if malfunction can not be removed by consumers themselves

Contact Us



MALAYSIA

Healthy World Lifestyle Sdn Bhd
No.22 Jalan Anggerik Mokara 31/47
Kota Kemuning, 40460 Shah Alam
Selangor Darul Ehsan, Malaysia
Tel : 1800 88 4688



HONGKONG

OGAWA Health Care
International (HK) Ltd
31/F, YHC Tower,
1 Sheung Yuet Road, Kowloon Bay,
KLN, Hong Kong S.A.R.
Tel : +85 2 3465 6262
Fax: +85 2 2410 0119



SINGAPORE

OGAWA Healthy-Care Pte Ltd
27 Tampines Industrial, Avenue 5
T5 @ Tampines
Singapore
Tel : +065-6846 8588
Fax: +065-6846 7677



VIETNAM

OGAWA Vietnam
2nd Floor , Petroland Tower ,
12 Tan Trao , Tan Phu Ward ,
District 7 , Ho Chi Minh
Tel : +84 8 3775 3093 / 3775 3094
Fax: +84 8 3775 3095



PHILIPPINES

O-Healthcare Solution Phil. Inc.
2303, 4th Floor, G&A Building,
Chino Roces Avenue Extension,
Makati City, Philippines
Tel: +632 889 0676 / +632 844 6851
Fax: +632 812 0421

目录

安全指引	14
保修权限	16
产品规格	16
产品的部件名称说明	17
使用操作说明	18
自动程序	19
使用姿态	21
产品的清洁与保养	22
常见问题及处理方法	23
联系我们	24

安全指引

请在使用产品之前仔细阅读此使用说明。

本产品仅使用电压为220-240V ~，频率为50/60Hz的电源。

i. 使用环境

- 为避免机器故障或者触电，请勿在潮湿或者多尘的环境中使用本产品。
- 请勿在浴室使用本产品。
- 请勿在床上使用本产品。
- 请将本产品放置在安全的位置，避免浸水。
- 产品放置位置请远离火源或热源，避免阳光直晒。
- 请在平坦地面使用本产品。

ii. 健康预防

此产品不适宜以下人士使用：

- 心脏病患者、恶性肿瘤患者、背痛、脊椎侧弯或异常、高血压、骨质疏松症等；
- 孕妇或月经期间；
- 身体不适者、发烧、身体有红肿或发炎、瘀青、过敏等；

iii. 以下情况下请勿使用本产品

- 连续使用本产品20分钟以上；
- 孩童或者行动不便者在无监护的情况下请勿使用本产品；
- 困乏或者醉酒时；
- 液体不慎溅落在产品上时；
- 使用中感觉任何不适或者疼痛；
- 使用中发现机器有任何故障或者异常；
- 打雷时

安全指引

iv. 警告

- 请将插头完整插入插座后方可使用。
- 没有使用时，请拔下电源插头。
- 请勿使用非OGAWA提供的配件。
- 如遇本产品电源线脱落或者表皮破损，切勿使用。请联系OGAWA客服中心进行检修。
- 勿将任何物体掉入或者插入产品机构间隙。
- 请勿扯动电源线来移动本产品。
- 请勿在运行时遮盖本产品。
- 请勿挤压或用尖锐物划伤本产品。
- 体重超过100公斤以上禁止使用。
- 本产品仅供个人使用，非用于商业用途。
- 本产品每次仅限一人使用。

v. 注意事项

- 肢体、感官或精神上有残障或缺乏相关知识经验的人必须在被告知相关风险以及正确使用方法，并有他人指导监督的情况下才可以使用该产品。
- 过度使用可能会导致产品过热而缩短产品寿命，连续使用时，须停用一段时间让产品冷却后再继续使用。
- 正在接受药物或医疗治疗者，请在使用前咨询医师。若您感觉不适，请勿使用本产品。
- 请勿拆卸本产品。触碰内部机件可能导致故障或触电意外。不可试图自行拆卸或者修理本产品，本产品仅由专业维修人员保养修理。
- 请勿用潮湿的手安插或拔出电源插头。
- 当电源插座有松动时，请勿使用。



保修权限

- 1. 本保修权限不可转让。
- 2. 本产品仅对购于OGAWA或其授权代理的产品有效。
- 3. 当要求保修服务时，请出示购买凭据（正本收据）。
- 4. 所有保修项目，包括但不限于这些健身及适销性的保修项目，其有效期限需视保修政策而定。
- 5. 因错误使用，处理不当以及未经授权私自改装或维修产品导致的损坏不属于保修范围。
- 6. 保修范围不包括不属于产品本身的附件。
- 7. 保修不适用于租赁、交易、商业、机构及其他非居民使用者。
- 8. 保修政策涉及的维修服务须经OGAWA授权，且只有获授权的技术人员才可进行维修。
- 9. 若需要更换的维修附件缺货，OGAWA有权采用替代配件进行维修。

产品规格


产品型号	OE 1031
额定电压	220-240V~
额定频率	50/60Hz
额定功率	200W
额定工作时间	10 minutes
产品重量	6.5KG
产品尺寸(长*宽*高)	367x282x127mm

产品的部件名称说明



使用操作说明

控制器的名称和功能



1. 电源按钮：打开和关闭。控制面板将亮起并将进入待机模式。默认程序为手动模式；屏幕显示“88”。

2. 速度 +/- 按钮：选择手动模式后，按按钮增加/减少振动速度，从 1-30 挡。

3. 速度快捷键：10 挡, 15 挡, 20 挡。

4. 自动程序：P1, P2, P3

5. 时间 +/- 按钮：选择手动模式后，按按钮增加/减少时间，从 1 到 10 分钟。

6. 启动/停止按钮：按下按钮启动，再次按下停止。



7. 模式按钮：按下按钮切换 3 个预设程序，P1-P3 和手动模式。

*预设无法调整速度和时间




备注:远程控制信号的有效接收范围为2.5米，控制器的红外发射指示器应与产品的红外接收窗口对齐。

自动程序



自动程序



使用姿态

A



B



C




D



E




F



G



H



I



产品的清洁与保养

清洁之前，请确认电源开关已经关闭，电源插头已经从插座拔出。脚踏皮和外壳可用含清洁剂的湿布彻底清洁，然后用干布小心擦除清洁剂。



不要用苯或稀释剂等腐蚀性洗涤剂清洗产品，因为它可能会导致产品的褪色、腐蚀或开裂。不要把水溅到产品上以避免短路。

控制器

由于湿布可能会引起控制器的故障，只能使用干布清洁。

储存维护

如果产品不使用时，确保电源开关关闭电源插头已经切断的情况下清除产品上的灰尘和污垢，保持清洁，远离高温或高湿度的地方。

保险丝的更换

将熔断器的方向转向箭头方向(逆时针)，用顺时针方向的螺丝把保险丝更换。

常见问题及处理方法

S/N	现象	可能原因	处理方法
1	无法开机	电源插头没有插好	插好电源插头
		产品电源开关未打开	打开电源开关指示灯亮
		瞬间电流过大，保险丝熔断	更换同型号的保险丝
2	机器工作中发出"咻咻"或"哒哒"的声音	这是机器运转时发出的声音，而非机械故障，属正常工作声音	无需处理
3	遥控器操作，机器无反应	遥控器电池没电	更换电池
		遥控器遥控距离太远	无线感应的有效范围为3M
4	无法正常操作	可能执行命令过于频繁	关闭电源，重新启动
5	显示“ER 02”	通讯故障	关闭电源3分钟后再开机，如不行请联系售后
6	显示“ER 01”	产品过载保护	关机30分钟后，再次使用

备注：当消费者无法自行排查原因时，请联系售后。

联络我们



马来西亚

Healthy World Lifestyle Sdn Bhd
No.22 Jalan Anggerik Mokara 31/47
Kota Kemuning, 40460 Shah Alam
Selangor Darul Ehsan, Malaysia
Tel: 1800 88 4688



香港

OGAWA Health Care
International (HK) Ltd
31/F, YHC Tower,
1 Sheung Yuet Road, Kowloon Bay,
KLN, Hong Kong S.A.R.
Tel : +85 2 3465 6262
Fax: +85 2 2410 0119



新加坡

OGAWA Healthy-Care Pte Ltd
27 Tampines Industrial, Avenue 5
T5 @ Tampines
Singapore
Tel : +065-6846 8588
Fax: +065-6846 7677



越南

OGAWA Vietnam
2nd Floor , Petroland Tower ,
12 Tan Trao , Tan Phu Ward ,
District 7 , Ho Chi Minh
Tel : +84 8 3775 3093 / 3775 3094
Fax: +84 8 3775 3095



菲律宾

O-Healthcare Solution Phil. Inc.
2303, 4th Floor, G&A Building,
Chino Roces Avenue Extension,
Makati City, Philippines
Tel: +632 889 0676 / +632 844 6851
Fax: +632 812 0421



Ogawa is grateful
to have received
these awards:



PUTRA
BRAND AWARDS 2010
SILVER
In association with
MALAYSIA'S MOST VALUABLE BRANDS

Ogawa Wellness Hubs: Malaysia.Singapore.Vietnam.Philippines.Indonesia.Myanmar.India.China
Hong Kong.Taiwan.Turkey.Pakistan.Mauritius.Russia.France.Canada
Thailand.Korea.Au stralia.Dubai.K.S.A.Nigeria.Egypt.Lebanon.Denmark
Sweden.U.S.A