
User Manual

WonderEase

Portable Neck and
Back Massager

OL0278

Please read this manual carefully before use and
keep it in a safe place for future inquiry.

IMPORTANT SAFETY INSTRUCTION

When using the product, all instructions and basic precautions should always be followed, including the following:

DANGER - To reduce the risk of electric shock:
Always unplug this product from the electrical outlet immediately after use and before cleaning

WARNING - To reduce the risk of burns, fire, electric shock, or injury to persons:

- 1) The product should never be left unattended when plugged in. Unplug from outlet when it is not use and before putting on or taking off parts.
- 2) Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- 3) Never operate the product with the air openings blocked. Keep the air openings free of lint, hair, and the like.
- 4) Do not drop or insert any object into any opening.
- 5) Do not use outdoors.
- 6) Do not operate where aerosol (spray) products are being used or where oxygen is administered.
- 7) Be sure to pull the plug out of the socket when you reach in. And don't pull it out when your hands are wet. That could cause electric shock and injury.
- 8) Do not carry this product by its supply cord or use cord as a handle.
- 9) This product contains batteries, it can only be replaced by technicians.
- 10) Remove the power plug from the socket immediately after the power failure. Otherwise, it may cause electric shock, short circuit and fire.
- 11) Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the appliance to a service center for.

-
- 12) Special attention to components most susceptible to wear an advice to replace defective components immediately and/or keep the equipment out of use until repair.
 - 13) Use this appliance only for its intended use as described in this manual.
 - 14) Keep the cord away from heated surfaces. If the fault is found, stop using immediately and cut off the power supply. Avoid accidents or injuries. please do not repair it by yourself, please contact the seller or manufacturer in time, avoid accidents or injuries.
 - 15) For charging use only protection class II power supply approved for use with household product, SELV output 5V, max. 2A
 - 16) The appliance must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
 - 17) This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
 - 18) Children should be supervised to ensure that they do not play with the appliance.
 - 19) This appliance contains batteries that are non-replaceable.
 - 20) The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.

Please check whether the product and accessories are intact and damaged before using.

*** IF YOU FEEL ABNORMAL DURING USE,
PLEASE STOP USING IMMEDIATELY AND
CONSULT YOUR DOCTOR. OTHERWISE
MAY IT BE HARMFUL TO HEALTH.**

1. People who are susceptible to electromagnetic malfunctions such as cardiac pacemakers and implantable medical electronic devices.
2. People who are supposed to lose their sense of temperature.
3. People with malignant tumors.
4. People with heart disease.
5. A person who is restless early in pregnancy or shortly after giving birth.
6. One who has an abnormal or curved spine.
7. A person with a disturbance of perception caused by a high degree of peripheral circulation in diabetes.
8. Someone with a scar on their skin.
9. People who need quiet and people who are in poor health.
10. People with osteoporosis, severe fractures, sprain, and muscle pain, etc., acute (head) pain.
11. People with 38 °C temperature (fever).
12. (acute inflammatory symptoms, tiredness, chills, changes in blood pressure, etc.) a time of intense, debilitating stress.
13. People with tendonitis.
14. Others treated at medical facilities.

The device is not to be used for commercial use or medical purposes. If you have health concerns, consult your doctor before using this product.

Only use indoors!

Do not use the device in damp rooms (e.g. when bathing or showering)

- Do not wash!
- Do not dry clean!
- After using this product for 10 minutes, rest for 30 minutes, then use again.


SAVE THESE INSTRUCTIONS

Correct Disposal of this product

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

NAME AND FUNCTION OF THE MACHINE PARTS


A.Massage Node


B.Power Switch

C.Buckle & Handle

D.Type-C charging port on the back

Type-C
Charging Cable

HOW TO OPERATE

Power/ Mode

1. Long press the button for 2 seconds to start machine. The default mode is kneading and the light at white.

Long press the button for 2 seconds again, to closed the machine.

Automatic shutdown: The machine stops working after 10 minutes of use.


2. When machine is started, the default mode is kneading mode. Press this button to change the mode at kneading mode → Shiatsu mode → kneading mode in turn.

Heat

When machine is started, the default mode is not heat function.

Press the button, to start the heat function and the light at red. Press again to closed the heat and light at white in turn.

MASSAGE POSITION

Cervical massage

Lower back relaxation

Relief leg pressure

The shoulder-back design completely frees hands

Change the previous handle design, you can massage while working

HOW TO CHARGE

1. The machine cannot work when it's charging.
2. The machine is charging, the light is red.
3. The battery with full charging, the light is white.
4. When the machine's battery power is less than 20%, the light at red flash, please charge it on time.

※ It is normal if the unit and charging plug being heat up slightly during charging.

Charging Time

The charging time is approximately 3.5 hours.

After Charging

After charging, unplug the power plug.

Attention

- * When using this product, be sure to use the company's certified switch power cable.
- * This product is only equipped with Type-C charging cable for product charging, you need to buy your own charging head. The charging head is recommended to use DC 5V-2A .
- * If the product is not used for a long time, please charge the product within half a year to activate the battery

HOW TO CLEAN

Cleaning

Before you clean, please make sure the power is switched off, and the power plug has been pulled out from the socket. Please use dry cloth to clean the product.


Maintenance

Keep the product clean and far from the wet environment.

Clean when you feel cleaning is necessary.

Cleaning the main unit

- Wipe away dirt with a cloth soaked in soapy water, wring out the cloth before wiping.
- Do not use alcohol, paint thinner, benzine, or nail polish remover. Using these may damage the appliance or cause cracks or discoloration in the components.

PRODUCT SPECIFICATIONS

Model No	OL0278
Product Name	WonderEase Massager
Power output	DC 5V --- 2A
Power consumption	18 W
Battery capacity	2000mAh

OPERATIONS FOR USE


Problem	Cause	Remedy
Unable to start the device.	The cord is not connected properly.	Please push it harder until you heard the sound of "click".
The appliance stops during use.	The appliance automatically stops after approximately 10 minutes.	This does not indicate a problem.
The main unit does not heat up.	The appliance may not work beyond the recommended temperature.	Use under the recommended temperature of 15 °C - 35 °C

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