
Product Manual

Please read all of the operating instructions carefully before use and keep in a safe place for future reference.

Operating Instructions:

1. Long press to power on
2. Long press to shut down
3. When turned on, short press the on/off button to reveal the ambient light
4. Short press the on/off button again to turn off the ambient light
5. Low battery prompt blue light keeps flashing (power too low automatic shutdown); charging status red light keeps flashing (full charge status long green light)

Product parameters	
Product Name: On-the-Go Mini Massager	Product Model: OZ3022
Product Size: 58×58×27mm	Battery Material: rechargeable lithium-ion battery
Main Material: ABS+silicone	Battery Capacity: 400mAh
Light Function: 265 color ambient light	Net Weight: 47g
Timing: 15 minutes	Rated Power: ≤1.5W
Charging Time: About 1 hour	Input Voltage: DC5V/0.35A
Use Time: about 6 hours	Mode Selection: 12 kinds of massage modes automatic cycle
Executive Standard: GB4706.10-2008	(Vibration, Tapping, Tuina, Kneading, Scraping, Acupuncture, Rolling, Squeezing, Wipe, Rubbing, Palm Rubbing, Finger vibration)
Patent No.: 202330278646.6	

Matters needing attention

⚠ Danger

- Please do not use with any of the following medical electronic devices:
 1. Implantable medical electronic devices in the conductor of cardiac pacemaker;
 2. Artificial heart and lung and other life-preserving medical electronic devices;
 3. Medical electronic devices such as electrocardiograph;* may cause the misoperation of these medical electronic devices, which may endanger life;
- 4. The following people are contraindicated: pregnant women, patients with abnormal heart or blood pressure

⚠ Attention

- Do not plug and remove the USB charging cable with wet hands. Failure to do so could result in an electric shock or accident.
- Do not knock or drop this product, repeated knocks may cause internal damage.
- Do not give it to the unconscious. Otherwise, it can cause accidents or physical discomfort.
- Do not use it on the same part for a long time (more than 30 minutes). Otherwise it may make the muscles of the massage part tired and cause physical discomfort.
- If you feel unwell due to this product, please stop using it immediately.

Troubleshooting

- If you encounter problems when using the product, please refer to the following tips and solutions.
- If the product is still in malfunction or unusable, please turn off the power.
- Please do not disassemble or repair the product by yourself. Damage to the product due to self-disassembly or repair will result in invalidation of the warranty. The Company will not be liable for any consequences arising from this situation.

Condition	Possible cause/handling methods
Unable to start	● Whether the internal rechargeable battery is low. ● Please confirm whether the start key of the mainframe product has been pressed.
Cannot charge or no charge light display	● Please check whether the USB cable and the host and power adapter are properly connected.
Display normal but no sensation	● Is the massager in direct contact with the skin? Adjust the position so it fits against the skin.
Tingling	● Q: Does the massage take too long? A: It is recommended to press for 15 minutes twice a day. ● Q: Is there any wound on the skin surface? A: Do not put electrodes on the wound for massage.
Stimulation is very weak	● Q: Does the massager fit snugly against the skin? A: Need to adjust the position so that it is close to the skin and adjust the appropriate strength. ● Q: Is the battery low? A: If the battery is low, please charge it.

Warranty Terms


- ★ The following cases are not covered under warranty:
 - Loss of the warranty card or expiration of the warranty period.
 - Damage caused by human factors, natural disasters, or unauthorized disassembly and repairs.
 - Malfunctions caused by improper use not following the user manual.

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