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OGAWA

by OGAWA

Curve EaseX
Extended
Massage Gun

OLO178

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5

User must read

Please read this manual carefully before using the product.

Symbol Description

To avoid personal harm or property damage, please be sure to observe the following:

· The contents desired to be observed are respectively explained by the following diagrams.

⊘ Indicates prohibition

① Indicates the need to follow suit.

Applicable crowd

①The following groups should be used with caution

· patients with heart disease

· malignant tumor patients

· Injured and not healed

· People with symptoms of fever

· People with muscle damage

· People during pregnancy or shortly after childbirth

· People with sequelae due to traffic accidents, etc.

· Need quiet, as well as poor physical condition of the crowd


· People with perception disorders caused by diabetes

· People with cutaneous sensory disorders, or who are presumed to have lost the sensation of temperature

1

OGAWA

by OGAWA

Curve EaseX
Extended
Massage Gun

OLO178

User's Guide

>>>> Please read the instructions carefully before use

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Method of use

Preparation before use

1. Please check the machine and accessories are complete.

2. Please remove the pendant and decoration when using to avoid affecting the massage part.

5

Use Operation

1./On/Off/:
Long press the switch button 2 seconds to boot: boot state, long press the button 2 seconds to fully charge the battery.

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Charging method

1. When using for the first time, please fully charge the battery. It takes about 3 hours to fully charge the battery.

2. need to use 5V=2A adapter for charging.

5V=2A adapter

Type-C charging port

6

Maintenance Maintenance

Cleaning method

①The silicone sleeve of the product cannot be removed.

1. The surface of the equipment can be cleaned with a slightly damp towel and then dried with a soft cloth.

2. Do not allow water to enter the inside of the product, please leave it to dry naturally.

① When this product is not used for a long time, please charge the product once every 1 month.

① Before maintenance, the charging cable must be pulled out, but do not use wet hands to pull out the charging cable.

Precautions for maintenance

① Be sure to turn off the machine when cleaning.

① Please properly handle the packaging bag of this product to keep it away from infants and young children.

① Maintenance, please do not use volatile oil, thinner, alcohol, strong alkaline detergent, bleach and other cleaning.

① Frequently check the fuselage for signs of wear or damage. If there are such signs, or the machine has been misused, you should contact customer service for maintenance before using it again.

① This product contains a lithium battery, please do not disassemble, modify, flood the product (controller) or throw it in the fire, otherwise it may cause damage to the product.

① Do not damage, modify, bend or pull the charging cable. In addition, do not place heavy objects on the charging cable or squeeze the charging cable.

7

① This product is suitable for cars, trains, ships, aircraft transportation. During transportation, the product should be fixed to prevent severe vibration, collision, rolling and falling, and put up, light and rain-proof.

Abandoned way

The lithium-ion battery used in this product is a valuable resource that can be recycled. Please take out the lithium-ion battery and hand it over to a professional institution that can recycle and regenerate the rechargeable battery.

Battery removal steps are as follows:
- Cut the connection between the handle and the silicone sleeve with a professional tool.
- Use a screwdriver to twist the handle screw;
- Twist the screws on the circuit board;
- Cut the battery lead wire with scissors;
- Remove the battery.

① Users cannot replace the batteries themselves.

The battery must be removed from the appliance before discarding the appliance.

The appliance must be de-energized when the battery is removed.

Batteries should be handled safely.

8

Fault diagnosis

This chapter lists the most common questions about this product. If you cannot find the answer to your question, please contact the service hotline or an authorized service center in your area.

Fault phenomenon

Possible causes

Processing method

Host cannot be charged

Internal rechargeable battery low

Host charging

Recharge the host

Host cannot be charged

Charging cable or adapter not connected properly

Reconnect

Suddenly stop at work

Internal rechargeable battery low

Host (safety protection) causes shutdown

Recharge the host

Suddenly stop at work

Host (safety protection) causes shutdown

Recharge the host

After a long period of non-use, be sure to confirm that the function is normal before you can use it again.

Whether the following faults:
· heating of charging cable or plug
· The charging cable has obvious scale or deformation, and has a burning smell.
· After the switch is turned on, it does not work properly
· Other abnormalities or malfunctions

If the above fault occurs, in order to prevent accidents, please immediately stop using, and cut off the power supply. Please send the maintenance point of our company.

9

Main technical parameters

Product Name

Curve EaseX Extended Massage Gun

Product Model

OLO178

Product size

368.6×210×48 mm

Product Weight

About 955g

Battery capacity

1500mAh

Operating voltage

7.4V

Working power

15W

Charging voltage

5V=2A

Maximum travel

10mm

Interface Type

Type-C

Material

ABS、Silicone

Product List

Host

x1

Massage head

x4

Charging line

x1

Instructions

x1

Certificate

x1

10

Environmental Instructions

1. This product meets the requirements of the management method of restricting the use of harmful substances in electrical and electronic products;

2. In order to protect the environment, after the product is scrapped, as a consumer, you have the responsibility to separate it from the domestic waste and send it to a qualified recycling station, which will carry out classification, disassembly, recycling and reuse in accordance with relevant national regulations.

Name and content table of hazardous substances in the product

Part name

(Pb)

(Hg)

(Cd)

(Cr(VI))

(PBBs)

(PBDEs)

Battery pack

○

○

○

○

○

○

Wire components

○

○

○

○

○

○

Plastic components

○

○

○

○

○

○

Charging cable

○

○

○

○

○

○

○: Represents that the content of this hazardous substance in all homogeneous materials of the component is below the limit requirements specified in GB/T 26572-2011.

×: It means that the content of the hazardous substance in at least one of the homogeneous materials of the component exceeds the limit requirements specified in GB/T 26572-2011.

10

Under the normal use conditions described in this specification, the environmental protection service period of the product is 10 years, which can be safely used within the 10-year environmental protection use period. After the environmental protection service period, it should enter the recycling system.

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