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
Important Safety Instructions

When using this appliance, basic precautions should always be followed, including the following:

Read all instructions before use.

1. Do not operate this product with a damaged USB charging cable or any damage to the product. To avoid electric shock, please do not disassemble and reassemble. If there is any damage, please bring it to our repair centre. Unauthorized repairs/modifications or products used for purposes other than the original design of the product will not be covered by the warranty.
2. Do not allow children to use this product alone.
3. Keep this product away from water, fire and harsh environments.
4. If the product falls into the water accidentally during charging, please unplug the USB charging cable immediately, and do not directly touch the product that has fallen into the water.
5. Do not use this product in humid, high temperature and high electromagnetic environment.
6. This product contains lithium batteries, please do not disassemble, modify or throw in the fire to avoid danger.
7. This product has heating function, please pay extra attention while using it.
8. Do not use this product continuously, otherwise, it may cause fatigue and discomfort on the massaged part.
9. This product is not suitable for any patients or people who are incapable of operation.
10. This product is a household electronic product, not as a medical device, and cannot be used for diagnosis, treatment and other commercial purposes.

Special statement:

1. Please keep this manual for future reference.
2. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safety.
3. Children should be supervised to ensure that they do not play with the appliance.
4. Do not immerse the product, the USB charging cable or the magnetic charging point of the product in any liquid. Otherwise, it may cause electric shock or fire. To clean this product, please follow the cleaning part of the instruction manual.
5. When the product is not in use, please turn it off.
6. Do not touch the product with wet hands.
7. Do not expose the product outdoor (rain or under the sun).
8. Before using this product, remove all packaging materials.
9. After unpacking, check for damage and there are no missing parts.

10. All packaging materials, such as plastic bags, polyfoam, tape, etc. should be placed in places where children are not exposed because improper use of these materials will lead to child suffocation and even death.
11. Use only the parts provided in the product and the packaging, otherwise, it may pose a hazard to the user and may damage the product.
12. Before using the product, please check whether the power and voltage are consistent with the technical specification of the product. Technical specification is shown on the gift box and this instruction manual. If you have any questions, please contact our customer service centre for assistance.
13. Do not allow the USB charging cable to touch a sharp edge or hot surface during operation.
14. This product is for home use only, not for any commercial purposes.
15. If the product needs to be repaired, client must clean the product before sending to our repair centre, otherwise, our company will charge extra cost for cleaning.

Warning:

Do not use this product if there is one of the following conditions:

- A. Pregnant women and those with heart disease and medical history.
- B. Those who have recovered from abdominal or lumbar surgery.
- C. Patients who are undergoing treatment and feel abnormalities.
- D. Patients with hypertension, arteriosclerosis or cerebral thrombosis.
- E. Patients with high fever, malignant tumor and osteoporosis.
- F. Patients with infectious diseases.
- G. Patients with skin sensation disorders or skin abnormalities.
- H. Those who are uncertain about their health.

Package Instruction


Before using the product, make sure the accessory is complete. In case of missing or damaged items, please contact Ogawa and replace them with new ones.

The list of accessories in the package is as follow:

1. BellaX Slimming & Heating Belt With Vibration Massager x 1
2. Magnetic charging cable x 1

Function icon:

- | | |
|---------------------------|----------------------------|
| 1. On/Off mode | 5. Vibration adjustment |
| 2. EMS mode | 6. Magnetic charging point |
| 3. EMS higher strength | 7. EMS lower strength |
| 4. Temperature adjustment | 8. Strap eject button |

Function operation instructions

1. On/Off Mode

Long press 3 seconds to turn On/Off the machine, short press once to switch the mode. This product has two modes; Warm mode and Slimming mode.

2. EMS mode key:

In the Slimming mode, short press once to switch the EMS mode. There are three modes: Hammer, Kneading and Dynamic mode.

3. Temperature adjustment key:

Short press once to switch the temperature. There are three temperature levels: low, medium, high temperature, and temperature shuts off (the temperature shuts off in Slimming mode only).

Note: The low temperature is approx. 39.5°C, the medium temperature is approx. 41°C, and the high temperature is approx. 43.5°C (the temperature difference between each level is about +/- 1 degree).

4. **Vibration adjustment key:**

In the Slimming mode, short press once to switch the motor vibration. There are levels of vibration frequency: slow, medium, fast.

5. **EMS strength plus and minus keys:**

In the Slimming mode, short press once to increase or decrease the intensity. There are ten levels of EMS intensity.

Note: The EMS pulse and vibration functions can only be turned on in Slimming mode; in Warm mode, only red light and heating functions are available.

Preparation before use

Please confirm whether the product is sufficiently charged. If the battery is too low, the machine function will turn off or cannot be turn on, please charge it before using. After the battery is fully charged, be sure to unplug the USB charging cable.

Recharge

1. Connect the USB charging cable to the magnetic charging point of the product.
2. The indicator light is red during charging, and blue after charging is complete.
3. During the charging process, the product does not support booting, the product starts the self-protection program and stops working.
4. The appliance adopts a built-in rechargeable battery, the charging time is 2-3 hours.
5. During use, when the battery is low, a voice will prompt to indicate that it is low battery. Please suspend use immediately and only use it after fully charged.

Multiple charging methods

Computer USB interface, mobile power supply, 5V /1A power adapter

Note: The default working time of power-on is 15 minutes, if you need to extend the use time, you need to turn on the host again.

Clean & Maintenance

- Before cleaning and maintenance, make sure to unplug the USB charging cable and the product has been completely cooled down.
- Please use a clean, dry soft cloth to wipe the main body. (**Note:** Please do not wipe the operation panel and the USB charging cable with a wet soft cloth.)
- Cannot use banana oil or thinner solution, gasoline, decontamination powder, hard brush, etc. to clean all parts of the appliance.
- Do not let water seep into the appliance. Do not put the electric parts in water.

Technical Specification

Product Name	BellaX Slimming & Heating Belt with Vibration Massager				
Model	Rated Voltage	Rated Current	Battery Capacity	Weight	Dimension
OM1130	5V 	1A	1200mAh	260g	225 x 115 x 45mm

Troubleshooting and maintenance

1. If a malfunction occurs during operation, please check the following. If the fault cannot be eliminated or the product still cannot operate normally, please turn it off and have it repaired by professional technicians. Do not disassemble or assemble the product by yourself.
2. If the product still cannot be resolved after troubleshooting, please call the customer service centre or repair centre.

Fault situation	Reason	Solution
Cannot boot	Low battery	Use after charging
Weak pulse intensity	Users have weak perception of intensity	Increase the strength before use
No pulse intensity output	The skin is dry or the skin is weak in pulse perception or the intensity is at level 1.	Adjust the pulse force to 2 to start / use after rubbing water or essential oil on the skin surface
"Beep" sound prompts when using EMS	Product electrodes are not in contact with the skin	Please fit the product tightly to the skin until the sound disappears

WARRANTY

Healthy World Lifestyle Sdn Bhd (hereinafter called "the Company") offers one year warranty service to the product (commencing from the date of purchase invoice).

The warranty includes Spare Parts & Labour Fee. Services beyond the said warranty period will be charged accordingly to the quotation from our Service Department.

1. Customer should follow the operating instruction in operating the product. The Company shall not be liable for any direct or indirect loss or damage to any person or property, howsoever caused arising from improper use, or any stoppage, breakdown or failure of the product.
2. The warranty service will not be offered if the model, serial number of the product or the invoice has been altered or effaced.
3. This warranty does not cover any damage or loss caused by man-made, misuse, neglect, unauthorized repair / modification or use of the product for any purpose other than the original designed.
4. For warranty service, customer(s) is required to send the respective OGAWA product(s) to any OGAWA outlet for warranty service. Shipping and Transportation fees will be borne by customer(s), if required.
5. All terms stated are subject to change without notice. In case of dispute, the Company reserves the right for final judgment and decision.

Contact Us

Healthy World Lifestyle Sdn Bhd (Malaysia)
No. 22 Jalan Anggerik Mokara 31/47
Kota Kemuning, 40460 Shah Alam
Selangor Darul Ehsan, Malaysia
Tel: 1800 88 4688