
ogawa
by OGAWA

((pulse))*plus*

User Manual 使用手册

Model: OE1032

CONTENTS

Safety Instructions	2
Warranty Claims	4
Specification	4
Name and Function of the Product	5
Remote Control Instructions	7
Program Mode	8
Use method	9
Storage and Maintenance	10
Surface cleaning	
Common Problems & Solutions	11
Contact us	12

Safety Instructions

Please read the manual carefully before using the appliance.
Power source for this appliance is 220–240V~ ;
alternating current frequency is 50/60Hz.

i. Usage environment

- To avoid malfunction or electric shock, do not use the appliance in a wet condition or dusty environment.
- Do not use the appliance in the bathroom.
- Do not use the appliance on the bed.
- Place the appliance in a safe location to avoid water immersion.
- Keep the appliance away from fire and heat, do not expose it to direct sunlight.
- Please use the appliance on a flat floor.

ii. Health precautions

If you are experiencing health complications or conditions as stated below, please consult your doctor before use:

- . Suffering from heart disease, malignancy, scoliosis, hypertension, osteoporosis, etc.;
- . Pregnancy or menstruation;
- . Fever, swelling, irritation, inflammation, congestion, allergy or any other health concerns

iii. When the product should not be used

- If the use exceeds 20 minutes at one time.
- This unit should not be used by children or invalids without adult supervision.
- Do not use this product when you are sleepy or drunk.
- If water is accidentally spilled onto the appliance.
- If the user feels abnormal or severe pain during the massage.
- If the user detects any malfunction or other abnormal conditions during operation.
- If there is lightning.

Safety Instructions

iv. Caution and warning

- Please completely insert the plug into outlet before use.
- Unplug the appliance from the outlet when not in use.
- Do not use any attachments that are not recommended by OGAWA.
- Do not operate this appliance if it has a damaged cord or plug. Please return the unit to OGAWA Service Center for examination and repairs.
- Do not drop or insert any object into any opening.
- Do not lift or carry this appliance by the power cord or any cord.
- Do not cover the appliance during operation.
- Do not crush the appliance and avoid any sharp folds.
- Anyone whose weight is more than 100KG should not use the appliance.
- This appliance is designed only for personal use, and not for professional purposes.
- The appliance can only be used by one person at a time.

v. Note

- The appliance is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge. They need to be supervised and instructed about its usage by a person responsible for their safety.
- Extensive usage could lead to excessive heating and a shorter product lifespan. Should this occur, discontinue the use and allow the unit to cool before operating again.
- If you are under medication or with a medical condition, please consult your doctor before use. Please do not use the product if you are not feeling well.
- Do not remove the cover of the appliance. Mishandling the internal components may result in malfunction or electric shock. Do not try to repair this appliance. The maintenance and repairs should only be made by technician.
- Do not connect or disconnect the power plug from the electrical outlet with wet hands.
- Do not use the appliance if the electrical outlet is loose.


Warranty Claims

- 1. Warranty is not transferable.
- 2. This warranty is effective only if the product is purchased from OGAWA or any authorized distributors.
- 3. Proof of purchase (original receipt) is required for all warranty repairs.
- 4. All implied warranties, including but not limited to these implied warranties of fitness and merchantability, are limited to the duration stated in the warranty policy.
- 5. Damage due to misuse, improper treatment and unauthorized modification and repair are not covered by this warranty.
- 6. Warranty does not cover accessories and attachments which do not belong to the product.
- 7. Warranty is not effective to rental, business, commercial, institutional or other nonresidential users.
- 8. All services covered by this warranty must be approved by OGAWA and repaired by authorized technicians only.
- 9. If replacement parts for defective materials are not available, OGAWA reserves the right to make substitution in lieu of repair or replacement.

Specification

Model No.	OE1032
Power source	220~240V~
Alternating current	50/60Hz
Power consumption	200W
Max Timer	10 minutes
Weight	7.5KG
Dimension(L*W*H)	437x334x138.6mm

Name and Function of the Product

Power Button

- 1.Switches between shutdown and standby modes.
- 2.When powered on or in operation, press the power button once to enter shutdown mode.
- 3.Press the power button again to switch to standby mode.
- 4.The button cycles between shutdown and standby mode in sequence. When the power switch is turned on, the device remains in standby mode by default.

Infrared Receiving Window

This window receives signals from the remote control for operation.

Name and Function of the Product

Mode Button

Switches between automatic and manual programs

When powered on, press the mode button to cycle through the following programs: P1, P2, P3, P4, P5, and Standby Mode.

Start/Stop Button

(Starts or stops the selected program in automatic or manual mode)

1. After selecting the desired program, press the Start/Stop button once to begin operation.

2. Press the Start/Stop button again to stop the current operation.

Time Increase/ Speed Increase Button

In Standby Mode: Press once to activate the time increase function in custom mode.

–The default maximum is 10 minutes.

–Each press increases the time by 1 minute, up to a maximum of 10 minutes (pressing further has no effect).

In Operation Mode: Press once to increase the speed level in custom mode.

–The default maximum is 99 levels.

–Each press increases the speed by 1 level, up to a maximum of 99 levels (pressing further has no effect).

In Automatic Program Mode: This button is inactive.

Time Decrease/ Speed Decrease Button

In Standby Mode: Press once to activate the time decrease function in custom mode.

The default minimum is 1 minute.

–Each press decreases the time by 1 minute, down to a minimum of 1 minute (pressing further has no effect).


In Operation Mode: Press once to decrease the speed level in custom mode.

–The default minimum is 1 level.

–Each press decreases the speed by 1 level, down to a minimum of 1 level (pressing further has no effect).

In Automatic Program Mode: * This button is inactive.

Remote Control Instructions

Mode Key

(Switches between automatic and manual programs)

- When powered on, press the mode key to cycle through the following programs in sequence: P1, P2, P3, P4, P5, and Standby Mode.

Power On Key

(Switches between shutdown and standby modes)

- When powered on or in operation, press the power button once to enter shutdown mode.
- Press the power button again to switch to standby mode.
- The button cycles between shutdown and standby mode in sequence.
- When the power switch is turned on, the device remains in standby mode by default.

Start/Stop Button

(Starts or stops the selected program in automatic or manual mode)

- After selecting the desired program, press the Start/Stop button once to begin operation.
- Press the Start/Stop button again to stop the current operation.

Time Decrease / Speed Decrease Button

- In Standby Mode: Press once to activate the time decrease function in custom mode.
- The default minimum is 1 minute.
- Each press decreases the time by 1 minute, down to a minimum of 1 minute (pressing further has no effect).
- In Operation Mode: Press once to decrease the speed level in custom mode.
- The default minimum is 1 level.
- Each press decreases the speed by 1 level, down to a minimum of 1 level (pressing further has no effect).
- In Automatic Program Mode: This button is inactive.

Time Increase / Speed Increase Button


- In Standby Mode: Press once to activate the time increase function in custom mode.
- The default maximum is 10 minutes.
- Each press increases the time by 1 minute, up to a maximum of 10 minutes (pressing further has no effect).
- In Operation Mode: Press once to increase the speed level in custom mode.
- The default maximum is 99 levels.
- Each press increases the speed by 1 level, up to a maximum of 99 levels (pressing further has no effect).
- In Automatic Program Mode: This button is inactive.

Note: The effective reception range of the remote control signal is 2.5 meters.


Program Mode

P2


circularly display up to 10 mins

P1


circularly display up to 10 mins

P3


circularly display up to 10 mins

P4

circularly display up to 10 mins

P5

circularly display up to 10 mins

Use method

Stand tiptoe in the middle of the massaging device. Straighten your back. Your abdominal muscles should be tight and you should feel changes in the tension in your calf muscles. You can also bend your knees 90 degrees. This stance will relax your waist and hips.

Standing on the massage device will align your feet and shoulders. Keep your back tight, your knees bent, and slowly tighten your leg muscles. To increase exercise intensity you should bend your knees as much as possible. You will begin to feel tension in your quads, hips, and back. This stance can exercise your thigh muscles, alleviate pain and soreness, promote blood circulation, and reduce fat deposits.

Kneel on the ground, place your hands on the product at shoulder width, and keep your body leaning forward. This stance will exercise your shoulder joints and arms.

Sit in front of the massager, with your back straight and your hands placed on either side. Place your calves on the massager, with your toes facing up. This stance should be used to completely relax while enjoying a powerful calf massage.

Sit on a chair in front of the massager. Place your feet on the device. The product's vibrations will help you relax the muscles in your back, buttocks, and thighs.

Storage and Maintenance

Surface cleaning

Before cleaning, ensure the power switch is turned off, and the attaching plug is pulled out of the socket.

The rubber pedal and covers can be cleaned with a wet cloth containing detergent, and the detergent should be removed with a Cloth

Notice

Do not clean the Product with corrosive detergent such as benzene Or thinner, as it may cause fading, corrosion or crack on the Product. Do not splash water into the Product to avoid short circuit.

Controller

As wet cloth may cause malfunction of the controller, only dry cloth Can be used.

Storage an Maintenance

Make sure that the power switch is turned off and the plug is removed From the power supply

Clean off the dust and dirt on the Product and keep it clean and away From places with high temperature or high humidity.

If the Product is not in use for along period, a dust proof cover (buy it Yourself) should be applied on the Product to prevent the entry of dust.

Common Problems & Solutions

S/N	Problems	Causes	Solutions
1	The Product can not be started	The plug is not in place	Insert the plug
		The power switch is off	Turn switch on, the indicator
		The fuse is burn	Replace with a new fuse of the same mode
2	Sound like phew or Da-da is heard during working	Normal sound generated by the Motor or from the abrasion Between mechanical units	No additional treatment is required
3	Operate the controller, but the Product no action	The electricity of batteries is worn out	Replace batteries
		The controller is too far away from the Product	The effective distance of infrared signal is 2.5M
4	Sudden stop	Default working time is up	Restart it after 10 minutes' interval
5	Can not operation	Too many commands sent to the Product and they can not be responded in time	Power off and then restart

Note: Please consult after-sales department, if malfunction can not be removed by consumers themselves

Contact Us

MALAYSIA
Healthy World Lifestyle Sdn Bhd
No.22 Jalan Anggerik Mokara 31/47
Kota Kemuning, 40460 Shah Alam
Selangor Darul Ehsan, Malaysia
Tel : 1800 88 4688

HONGKONG
OGAWA Health Care
International (HK) Ltd
31/F, YHC Tower,
1 Sheung Yuet Road, Kowloon Bay,
KLN, Hong Kong S.A.R.
Tel : +85 2 3465 6262
Fax: +85 2 2410 0119

SINGAPORE
OGAWA Healthy-Care Pte Ltd
27 Tampines Industrial, Avenue 5
T5 @ Tampines
Singapore
Tel : +065-6846 8588
Fax: +065-6846 7677

VIETNAM
OGAWA Vietnam
2nd Floor , Petroland Tower ,
12 Tan Trao , Tan Phu Ward ,
District 7 , Ho Chi Minh
Tel : +84 8 3775 3093 / 3775 3094
Fax: +84 8 3775 3095

PHILIPPINES
O-Healthcare Solution Phil. Inc.
2303, 4th Floor, G&A Building,
Chino Roces Avenue Extension,
Makati City, Philippines
Tel: +632 889 0676 / +632 844 6851
Fax: +632 812 0421

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Thailand.Korea.Au stralia.Dubai.K.S.A.Nigeria.Egypt.Lebanon.Denmark
Sweden.U.S.A